

MODUS[®] ESWT



*Superior Shockwave
Theraphy*

İnceler Medikal
www.incelermedikal.com
[CE 1984, ISO 13485]

*Where technology
and efficiency is delivered
with **style!***



MODUS is the most reliable,economical and effective ESWT System which you can reach.

- 1.MODUS has a very easy to use , portable, lightweight and very robust device.
2. Except compressed air,don't need any supplies to create cost.
3. MODUS hand applicator (Hand Paice) can be optionally 3,000,000 shots kit.
4. MODUS portable structure , can be moved easily if desired separated from the carrier table.
5. There is a suitable cubby for hand piece and a drawer for accessories on the side of MODUS.
6. Modus contains 7(seven) different kind of replaceable tip cap.(Transmitter- Applicator)
7. When the desired , treatment can be stopped and started with the help of trigger on hand applicator. (OPTIONAL)
8. There are expressions in MODUS software for all treatments, how and to which of the applications shold be done .
9. MODUS can be set automaticallly to stop the treatment when reaches the total number of shots for a session.
10. MODUS's pressure,frequency and all parametters such as the number of shots shock can easily readable on LCD screen.

MODUS-ESWT SET

Modus Main Unit	1 Unit
Handpiece	1 Unit
Applicators	7 Units

Ø 6 mm	P.Depth 0-45 mm
Ø 10 mm	P.Depth 0-45 mm
Ø 15 mm-Radial	P.Depth 0-45 mm
Ø 15 mm-Focus	P.Depth 0-45 mm
Ø 15 mm-Trigger	P.Depth 0-60 mm
Ø 20 mm	P.Depth 0-50 mm
Ø 35 mm	P.Depth 0-50 mm

Trolley	1 Unit
Compressor	1 Unit
Gel	1 Unit
Foot Pedal	1 Unit
Portable Carrying Bag	1 Unit

Major Therapuetic Areas

Plantar Fasciitis or Fasciosis
(Strained Arch)
Achilles Tendinitis or Tendinosis
Calcific Tendinitis or Tendinosis
Connective Tissue Pain
and degeneration
Muscle Pain and Injuries
Joint Injuries
Morton's Neuromas
Stress Fractures
Avascular Necrosis
(A dead portion of bone)
Slow-healing bone (Delayed unions)
Non-healing bone (Non-unions)
Rotator cuff tendinopathy
Epicondylopathy
Tibialis anterior syndrome

MODUS[®] ESWT



■ MODUS PROVIDES EFFECTIVE RHYTHMS AND FUNCTIONAL USE.

■ PORTABLE, PRACTICAL AND EFFECTIVE.

modus ESWT[®]

■ DOES NOT REQUIRE MEDICATION, ANESTHESIA OR SURGICAL INCISION.

**Quick results,
High success!**

Modus ESWT designed and manufactured by shockwave treatment expert Inceler Medical also provides easy treatment options with a wide range of applicators.



01

6 mm, Radial/0 - 45 mm
Acupuncture with Shockwave Therapy



02

10 mm, Radial /0 - 45 mm
Radial Shockwave Therapy, Tendinopathies



03

15 mm, Radial /0 - 45 mm
Radial Shockwave Therapy, Tendinopathies



04

15 mm, Trigger /0 - 60 mm
Chronicle Diseases, Deep Placed Diseases,
High Energy Needs



05

15 mm, Focus /0 - 50 mm
Head and Neck, Fascial Muscles,
Close to Surface Pains



06

20 mm, Radial/0 - 50 mm
Muscle Connective Tissue, Myofasciale Trigger Theraphy,
Trigger Diagnose, Skin Elasticity



07

35 mm, Radial/0 - 50 mm
Muscle Connective Tissue, Myofasciale Trigger Theraphy,
Trigger Diagnose, Skin Elasticity

MODUS ESWT TEKNİK BİLGİLER (Technical Specifications)

Manufacturer:	İnceler Medikal Sağlık Hiz. San. Tic. Ltd. Şti.
Model:	Modus®
Power Supply:	100 – 240 VAC, 50/60 Hz, 50 W
Fuses:	2 x 1 A, 250 VAC
Compressed Air Supply:	1 – 5 Bar (15 – 70 psi), 20 NI/min @ 5 bar
Frequency:	Up to 22 Hz
Operating Mode:	Single, Continuous, Burst
Display:	LCD
Classification:	EN-60601-1 Class I Applied part BF IP 20 (Foot Pedal: IP x 5) MDD 93/42 CEE Class IIb

Treatment Start / Stop: Main unit button, handle button and optional foot pedal

Fuse: O5 X 20 mm T4A 250 V

Air Pressure Support: 3.5 a 6.5 bar (3500 a hPa)

Body: Dimensions: Height: 120 mm
Width: 320 mm
Depth: 390 mm

Operating Conditions: Temperature: +10 to +30 C
Humidity: 30-75%
Atmospheric pressure: 700 - 1060 hPa

Storage and Transporting Conditions: -10 C to + 40 C
Humidity: 10% to 95%
Pressure: 500 hPa to 1060 hPa

Memory: 3, user defined

Custom treatment protocol selections: YES

[CE 1984, ISO 13485]

شرکت تین طب
021-88812365
www.teenteb.com/fa

